



# Patient Care Simulator "Cherry"




*Learning has never been more lifelike*

[limbsandthings.com](http://limbsandthings.com)

Introducing the new Patient Care Simulator, a basic full body manikin for training day to day patient care.

Designed for introductory nursing, life-like materials plus the ability to sit upright without support make this manikin ideal for realistic learning.



## Skills Gained

- Body positioning/patient handling/passive exercise
- Bed bathing and change of clothes
- Oral and denture care
- Nasogastric tube placement
- Oxygen/Nasal cannula placement
- Tracheotomy care and management
- Stomach pumping
- TPN care
- PEG care
- Stoma care
- Pressure ulcer care
- Subcutaneous injection
- Intramuscular injection
- IV insertion and care
- Male and female urinary catheterisation
- Enema
- Suppository insertion
- Perineal care



- BTEC and A Level, Health & Social Care courses within Further Education Colleges




- Male and female genitalia



## Product Qualities



- Seamless arms and legs enable realistic bed and foot bathing
- True-to-life joint articulation for patient handling and passive exercise
- Nasal and oral catheter can be inserted into the stomach




- TPN and IV routes to simulate patients with catheter lines
- Pressure sores of different depths located at the scapula, greater trochanter, sacral bone and heel
- Subcutaneous and intramuscular injection pads accept liquid and can be squeezed out after use



- Easily interchangeable male and female genitalia
- Manikin can sit upright without support
- Removable dentures



- Recommended tracheotomy catheter size ID8mm, OD11mm
- The PICC opening is 2cm depth
- Skin washable with soap and water



Product No: **KKMW22**

## You'll Receive

- 1 Cherry Patient Care manikin
- 1 Male genitalia
- 1 Female genitalia
- 1 Femoral IM injection pad
- 1 Upper arm IM injection pad
- 1 Wound care pad
- 1 Bottle
- 1 Syringe

Packed weight: 18.5 kg



Female and male genitalia



Syringe, bottle and IV injection pads



Wound care pad

## + Consumables + Replacements

- KK11401-010** Insulin Injection Pads (x6)
- KK11401-020** Gluteal IM pads (x4)
- KK11401-030** Female Genitalia Unit
- KK11401-040** Male Genitalia Unit
- KK11388-200** Injection Pads for Median Antebrachial Vein (x2)
- KK11401-050** Shoulder IM Pads (x2)
- KK11401-060** Thigh IM Pads (x2)
- KK11401-070** Wound Care Pad



*For demonstration or to order*

Limbs & Things Ltd. Sussex Street, St Phillips, Bristol, BS2 0RA, UK

T: 00 44 (0)117 311 0500 F: 00 44 (0)117 311 0501 E: sales@limbsandthings.com

**limbsandthings.com**