





PARACENTESIS TRAINER

This trainer fulfills the fundamental component of core medical training in diagnostic and therapeutic techniques of Paracentesis, for either palpation or ultrasound guided techniques.

Trainees are able to identify the echogenic anatomy enabling them to learn how to safely insert a needle or catheter into the peritoneal cavity.

60100 LIGHT

60111 DARK



Learning has never been more lifelike

limbsandthings.com

SKILLS

- Familiarity with the abdominal regions and underlying anatomy
- Palpation of anatomical landmarks
- Identification of excess fluid
- Through ultrasound, identification of safe needle/catheter insertion sites
- Insertion of needle into peritoneal cavity for therapeutic or diagnostic purposes
- Professional-to-patient communication

WHAT'S INCLUDED

60100

60111

Paracentesis Skin

Paracentesis Pads (x2)

Paracentesis Liver

Paracentesis Spleen

Paracentesis Bowel Insert

Paracentesis Floating Bowel

Paracentesis Clamping Frame

Paracentesis Large Volume Refill System

Base unit

Carry case



QUALITIES

ANATOMY

- Torso featuring bony landmarks (pubis symphysis, iliac crest) and umbilicus
- Internal anatomy includes:
 - Liver
 - Bowel
 - Spleen
 - Floating Bowel

REALISM

- Realistic tissue and needle response

KEY FEATURES

- Palpation or ultrasound techniques can be practised (side by side)
- Internal echogenic anatomy to allow recognition of landmarks under ultrasound
- Two 3.5 litre chambers can be filled with water for practise drainage
- Cost effective consumable pads

VERSATILITY

- Self-sealing pads, withstanding up to 200 needle or up to 100 rocket catheter insertions
- Ability to insert and remove a drain
- Positions both supine and lying on side

SAFETY

- Latex free

CLEANING & CARE

- Skin surface is washable using soap and water

CURRICULUM REFERENCES

UK

CORE MEDICAL TRAINEES:

Core Medical Training (CMT) ARCP Decision Aid, Joint Royal College of Physicians Training Board, 2015, p. 4

GASTROENTEROLOGY SPECIALIST TRAINEES:

Specialty Training Curriculum for Gastroenterology, Joint Royal College of Physicians Training Board, 2007, p. 38

ONCOLOGY SPECIALIST TRAINEES:

Specialty Training Curriculum for Medical Oncology, Joint Royal College of Physicians Training Board, 2007, p. 8

ACUTE CARE COMMON STEM TRAINEES:

Acute Care Common Stem Core Training Programme, Royal College of Anaesthetists, 2012, p. 87

ADVANCED NURSE PRACTITIONER TRAINEES:

Curriculum and Assessment Systems for Training in Emergency Medicine, Royal College of Emergency Medicine, 2015, p. 127

USA

EMERGENCY MEDICINE RESIDENTS:

University of Rochester Medical Center Emergency Medicine Residency Curriculum for Procedural Competency, 2011, p.4

INTERNAL MEDICINE RESIDENTS:

American Board of Internal Medicine Policies and Procedures, July 2017, p. 4

FAMILY MEDICINE RESIDENCY:

Procedures for Primary Care, p. 1380-1384.

SURGERY:

ACGME: Surgery FAQ, p. 5.

DIAGNOSTIC RADIOLOGY:

ACGME, Review Committee for Radiology, p. 1.

INTERVENTIONAL RADIOLOGY:

ACGME, Program Requirements in Interventional Radiology, p. 9.

International V.4 Mar 2021



REQUEST A QUOTE



REQUEST A DEMO



CUSTOMER SUPPORT

limbsandthings.com