

# BREAST EXAMINATION TRAINERS

Our Breast Examination Trainer provides a highly realistic learning platform for acquiring the skills required to perform Clinical Breast Examination (CBE).

Featuring 6 readily interchangeable and multi-positional pathologies, providing healthcare professionals with the tools to identify various complications and pathologies, including carcinomas, cysts, fibrocystic disease and fibroadenoma.

40201

41201

ADVANCED



This product was developed in close collaboration with leading clinicians:

- Mr Peter Donnelly, Torbay Hospital, UK
- Clinical Skills Resources Centre, University of Liverpool, UK
- Clinical Practice Centre, St James's University Hospital, Leeds, UK

40200

41200

STANDARD



40202

41202

BASIC



Both *simulated patient* and *bench top* training can be used in alignment with healthcare professionals promoting best practice techniques to trainees and patients, or for any undergraduate programs running **OSCEs**.

*Learning has never been more lifelike*

[limbsandthings.com](http://limbsandthings.com)



Soft tissue breasts look and feel realistic  
 Pathologies fit securely into supports  
 Clavicular and axilla pads for accurate lymph node placement  
 Breasts worn by a simulated patient  
 Self Breast Examination training

## QUALITIES

### ANATOMY

- Realistic soft tissue breast anatomy
- Pathologies supplied: carcinomas: 2cm, 3cm, 5cm, cyst, fibrocystic disease, fibroadenoma

### REALISM

- Soft tissue breasts look and feel realistic
- Clavicular and axilla pads for accurate lymph node placement

### KEY FEATURES

- Pathologies can be placed in various predetermined location points and are easily changeable
- Dual purpose product: bench top and Simulated Patient
- Hard torso available for bench top use

### VERSATILITY

- Comfortable to wear for long periods of time during OSCEs and assessments

### SAFETY

- Latex free

### CLEANING & CARE

- Skin surface is washable using soap and water

## SKILLS

- Clinical Breast Examination (CBE)
- Self Breast Examination (SBE)
- Identification of anatomical landmarks
- Identification of lymph nodes (axillary, supra & infraclavicular)
- Location and diagnosis of pathologies
- Professional-to-patient communication

## WHAT'S INCLUDED

	40201	40200	40202
	41201	41200	41202
Wearable Examination Breasts	✓	✓	✓
Pathologies	✓	✓	✓
Breast Examination Inserts	✓	✓	✓
Breast Pathologies Supports	✓	✓	✓
Breast Back Plates	✓	✓	✓
Breast Examination Torso	✓	✓	×
Lymph Node Pads and Lymph Nodes	✓	×	×

## OPTIONAL EXTRAS

40209 Medium Carry Case

### Patents granted

• US6336812 • EP099027 • DE69804793 • CA 2293585 • ZL98806331 • HK1023832 • 755575 Australia

USA V.2 Mar 2021



REQUEST A QUOTE



REQUEST A DEMO



BUY ONLINE



CUSTOMER SUPPORT

limbsandthings.com